

Digital bed and chair mats



1 About the Cassius Digital bed and chair mats

- The mats are placed on your bed or chair under a sheet or cushion.
- They have pressure sensors, and show when you sit or get up from your bed or chair.
- They also can show when you sit down or go to bed.
- Your carers can look at this information on their phone or computers and will get an alert if you get up or have not returned to your chair or bed.
- The mats can also be used on the floor or under a door mat to show your movement.
- As these mats are digitally connected so you need the internet for them to work.
- The mats are easy to wipe clean.



2 How does it work?

- The mats send an alert to your carer's phone by text message or to the portal if you get off your bed or chair. You can set a delay on the alerts to any time, for example 10 seconds or five minutes.
- Your carers can see information about your movements on the Alcove Online Portal and the Alcove Connect App. It will show them how long you are on/off the bed or chair and how long you spend in it.
- This will help your carers build a bigger picture in terms of times of need or to minimise risk.

3 Things to know

Internet

The devices need internet to be able to work and your carers will need access to the Internet or a smartphone and to connect to Alcove Online Portal / Alcove Connect App.

Battery

The sensor mats have a LS14250 battery, which has been designed to last between three and 15 years. It needs to be replaced by the user. If the battery is running low, you or your carer will get a text message notification from the Alcove system.

Charging

The mats are wirelessly linked to the black box, which needs to be connected to the power wall socket at all times. If it has been turned off the system will send a notification.

4 How is my data used for this product?

This device will capture data which is used to support best outcomes, you can find out more on our privacy statement - https://www.suffolk.gov.uk/about/privacy-notice/adult-and-community-services-privacy-notice/