



## The Cassius Falls Wearable

BUILT IN  
**FALLS**  
DETECTION

### 1 About the Falls Wearable

- The Cassius Falls Wearable is an attractive tear-drop shaped device which you can wear around your neck as a necklace or attached to your belt or keys.
- It will allow you to call for help if you need it – either at home or when out and about - and will give you the confidence to carry on with your day-to-day life.
- It also has a falls detector installed in it and will automatically raise an alarm should you have a fall.
- It is linked to a monitoring centre\* who will make sure you get the help and support you need, should you need it.
- It does not need internet to work but has an inbuilt SIM card.





## Overview



## 2 How does it work?

You can call for help by pressing and holding the two buttons together that are on the either side of the wearable in order to raise the alarm. A short pre-alarm sequence will start, and it will flash red.

An operator at the monitoring centre will then call your device and speak to you through the speaker and see how they can get the support you need.

Depending on how you fall, the device should detect this and automatically contact the monitoring centre without the need for you to press any buttons. If it doesn't do this, and you're able to, you will need to press the buttons together either side of the wearable to sound the alarm.

In the event of an emergency, your location can also be sent to the monitoring centre so support can be arranged – this is especially helpful if you do not know where you are.

## 3 Things to know

### Battery

The device includes a rechargeable battery and can be charged using the charger provided. The battery will normally last one to two weeks. The device will let you know by voice and will light up when it has 20 per cent charge left.

You can then place the device on the charger provided, which you can plug into the wall. The battery takes between three - five hours to recharge once placed on the charger.

### Waterproof

You can wear the device while washing in the shower or bath for a maximum of 30 minutes at one metre depth. It also works if it gets wet in the rain.

### Pacemaker

If you have a pacemaker then you should be aware that the device cannot be worn or come in contact within an area of 25cm from the pacemaker. If you are unsure then please check before use.

## 4 How is my data used for this product?

This device will capture data which is used to support best outcomes, you can find out more on our privacy statement - <https://www.suffolk.gov.uk/about/privacy-notice/adult-and-community-services-privacy-notice/>

\*The Monitoring Centre referred to in this sheet is managed by **Medequip**.